

SpectraPro v3.0x Software

Changes Since Version v3.3.5

1. Removed COM interface to enhance Vista compatibility.
2. Modified the internal v6 board setup motorized slit setup to include the new motorized slit.

Changes Since Version v3.3.4

1. Updated to bring the dll's in sync with updating SpectraSense from v4.4.1 to v4.4.2.

Changes Since Version v3.3.3

1. Updated to bring the dll's in sync with updating SpectraSense from v4.4.0 to v4.4.1.

Changes Since Version v3.3.2

1. Updated to bring the dll's in sync with updating SpectraSense from v4.3.9 to v4.4.0.

Changes Since Version v3.3.1

1. Updated to bring the dll's in sync with updating SpectraSense from v4.3.8 to v4.3.9.
2. Updated the SDK documentation (ARC_FilterWheel_dll.pdf, ARC_SpectraPro_dll.pdf, and ARC_Instrument_dll.pdf) to reflect additions to the dll.
3. Updated the SDK C header files (ARC_FilterWheel_dll.h, ARC_SpectraPro_dll.h, ARC_Instrument_dll.h) to reflect additions to the dll.
4. Updated the SDK Delphi Units (ARC_FilterWheelU.pas, ARC_SpectraProU.pas, ARC_InstrumentU.pas) to reflect additions to the dll.
5. The labview part of the SDK has been superceded by a new labview tool kit and is no longer on the CD.

Changes Since Version v3.3.0

1. Update to bring dll's in sync with updating SpectraSense from v4.3.7 to v4.3.8.
2. Added support for setting up the new Insight Spectrograph.

Changes Since Version v3.29

1. Added support so Mono Control Software could be launched from within SpectraSense.
2. Added support for Monochromators with custom Acton Sample Chambers (Wheel and XY Table Based). Note, this is not generic support for all sample chambers, but only custom motorized units from Acton Research Corp. that are attached to an SD-3.
3. Offset was not being correctly calculated for old model SP-150's. The issue has been resolved.

Changes Since Version v3.28

1. Fixed an issue with parsing the last grating in a grating list that prevented the user from using the last grating in the grating list on a SP-275. Issue may have been generic to the last grating position on all SP-xxx instruments.
2. Widened the instrument select dialog on the front page so the user can see the attached instruments full serial number.

Changes Since Version v3.27

1. Updated the dll's to keep in sync with SpectraSense V4.32

Changes Since Version v3.26

1. Modified the internal v6 board setup motor assign bitmap to show which side of the board the entrance slit is on.
2. Modified the internal v6 board setup hardware setup so that when a linear motor is selected mper is changed from 32 to 48.
3. Modified the internal v6 board setup assign motors functions so that the default min/max filter are 1 and 6.
4. Modified the internal v6 board setup default mper value for INS-300s to 48 from 32.

Changes Since Version v3.25

1. Added SpectraHub support
2. Synchronized the dll's with the current SpectraSense v4.30 release.
3. Fixed a bug in the internal v6 production version that did not correctly setup the motor app stop
4. Added support in the internal v6 production version to query and modify the board mper parameter.

Changes Since Version v3.24

1. Fixed a bug with it displaying the wrong current grating when one of the slots on a previous turret was blank.
2. Installed a work around for the early v6 instrument init-grating, turret bug.

Changes Since Version v3.23

1. Fixed a bug in the new internal v6 board setup of the side exit slit. It was assigning the wrong exit slit value.
2. Fixed a bug in the new internal v6 board setup display of the exit diverter assignment. It was displaying the wrong position for the diverter.
3. Added additional code so mono control software could be better used internally for production board setup.

4. Adjusted the decimal separator property so the software works better with European flavors of Windows.
5. Adjusted the time between the issuing of a ?mono-model and a who on an NCL to work around an issue older NCL's have with these two commands.

Changes Since Version v3.22

1. Updated the CD so it can install RS-232 Filter Control software.
2. Updated ARC_Instrument.dll to v0.0.7, Changes Included
 - a. Properly Jogs at Max Scan rate for all SpectraPro models.
 - b. Internal function to set v6 motor speed corrected.

Changes Since Version v3.21

1. Added ARC_Instrument.dll version info to the about screen.
2. Modified internal version to setup v6 board focal length, detector angle, and nm/rev.
3. Modified internal version to setup v6 board motorized indexable slit steps>microns, microns>steps, and home pointer.
4. Update ARC_Instrument.dll to v0.0.6, Changes include:
 - a. Added function to read an Instrument's (NCL, Mono, Filter, ect.) physical port number (COM1, COM2, ect.).
 - b. Added new function read nm per rev ration in linear instruments.
 - c. Added new internal v6 board setup function for setting nm per rev.
 - d. Added new Internal v6 function to read motorized slit clamp.
 - e. Added new Internal v6 function to read motorized slit freq.
 - f. Added new Internal v6 function to read motorized steps>microns, microns>steps, and home pointer.

Changes Since Version v3.20

1. VC++ and LabView examples added to the dll SDK.
2. Added support for Holographic Visible and UV gratings in the Blaze String.
3. Cosmetic touchup of how the grating list is displayed.
4. Added the COM Interface Active-X back to the Install, was omitted in v3.20 installer.
5. Added ability to restrict access to users changing grating parameters by creating a file called 'Restrict.File' and placing it in the "C:\Program Files\SpectraPro\bin" directory.
6. Full support for v6 controllers.
7. Open production factory setup parameters issues in the internal version for production of v6 units :
 - a. Read/Write nm per rev for linear systems
 - b. Read Slit Clamp Value
 - c. Read SFreq Value
 - d. Read Microns > Steps Value
 - e. Read Steps > Microns Value
 - f. Read Slit Init-Home Value

Changes Since Version v3.10

1. Now uses the correct grating offset when installing v6 mono controller gratings.-
2. Conversion to ARC_Instrument.dll, V3.20 and higher

Changes Since Version v3.09

1. Added support for the new v6 mono controller board.
2. Added support for filter wheels.
3. Modified the code that identifies active serial ports on a computer.

Changes Since Version v3.08

1. Corrected the software not finding valid AM/VM-510 monochromators.
2. Added additional manual support to the software.
3. Added support for a custom SP-275 with an SD-2 controller and custom slit assembly.
4. Modified how the init-wavelength and init-srate are handled.
5. Modified the internal behavior of the Splash Screen, it interfered with operations under Windows XP.
6. Upgraded the compiler from Delphi 6.0 to Delphi 6.0 Service Pack 2.

Changes Since Version v3.07

1. Tweaked the scan rate code. This way it takes into account the variable speed control feature of the motor in x00i monochromator.

Changes Since Version v3.06

1. Tweaked the G-Adjust and Offset code for SD-2 linear motors.
2. Clamped the wavelength limits in SD-2 controllers to reflect the stored firmware values.
3. Corrected the display of the initial scan rate being off by a factor of 10 when used with an SD-2.
4. Corrected the display of the initial wavelength being off for old model SD-748 controllers.

Changes Since Version v3.05

1. Made a modification that allows the software to work computer language settings that use a ‘,’ instead of a ‘.’ decimal separator.

Changes Since Version v3.04

1. Adjusted the timeout for motorized diverter mirrors in the slave unit of a double monochromator. The timeout was too short.
2. Added support for the X-Ray Grazing Incidence GIMS 551.5 monochromator.

Changes Since Version v3.03

1. Fixed an error involving the displaying of motorized slits when working with a double.
2. Updated a bug in the Wavelength calculators that showed up when using lasers of wavelength less than 100 nm.

Changes Since Version v3.02

1. Added support for the new Swift-X based prompts.
2. Added support for 1 and 2 gratings per turret in the x00i devices.
3. Modified the communications error handling for more robust communications.
4. Added the ability to abort a scan to keep it from repeating.
5. Added support for Doubles with motorized slits and diverter mirrors.
6. Fixed the slit not homing when the initialize button is pressed.

Changes Since Version v3.01

1. Added Holographic Blaze to the Install Gratings Page.
2. Added Mirror to the Install Gratings Page.
3. Added initial wavelength and initial scan support for the new swift-X firmware.
4. Corrected a problem with initial wavelength and initial scan rate for SD-2 based controllers.

Changes Since Version v3.00

1. Added a COM interface to the SpectraPro software.
2. Made changes necessary to support the COM interface.
3. Modified the WMCLOSE event, so it doesn't generate an error if nothing is attached.
4. Changed from the Roper logo back to the Acton Research Corp. logo.
5. Corrected the wording on the G-Adjust buttons.
6. Corrected a bug in which the job speed could be displayed at 10,000x what it actually was.
7. Added a Mono-Restore-Factory-Settings button.
8. Removed a blank restore defaults button on the Set Defaults page.
9. On the Set Defaults page, changed 'Initial Scan Rate' to 'Initial Scan Speed'.
10. On Set Defaults page, the 'Set New Values' now gives the user the option of canceling the operation.